Contec ProChlor Sporizid

pure¹¹-Nr.: 09111, Hersteller: CONTEC

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Quelle: <https://www.pure11.de/contec-prochlor-sporizid>

Contec ProChlor

Stabilised Hypochlorous acid blended with purified water

Contec ProChlor is a revolutionary new sporicide achieving a 100% kill against spores in under 1 min. A blend of hypochlorous acid in purified water, Contec ProChlor is provided ready-to-use and is efficacious against bacteria, fungi, moulds, yeasts, viruses and spores.

Contec ProChlor is 0.2 micron filtered, and filled under unidirectional air flow in a Grade B or C cleanroom. The sterile product is filled into pre-irradiated components. Provided double bagged, the product is designed for ease of entry into pharmaceutical cleanrooms.

Available as either a 1L trigger spray for small areas or in 0.5L, 1L and 5L capped containers for mop head saturation or when larger volumes are required.

FEATURES	BENEFITS
Log 6 kill against <i>Bacillus subtilis</i> in 60 secs	<ul style="list-style-type: none"> • Fast acting so saves time spent on biodecontamination • Fast enough to be used for transfer disinfection
Filtered to 0.2 microns in a Grade B/C cleanroom	<ul style="list-style-type: none"> • Ensures the product is free from contamination and particulates
Sterile version available	<ul style="list-style-type: none"> • Suitable for Grade A and B cleanrooms
Contains no quaternary ammonium or surfactant	<ul style="list-style-type: none"> • Very low residue, saving time on residue removal
No hazard classification	<ul style="list-style-type: none"> • Good operator acceptability as no strong odour • Only basic PPE required and no special disposal is required
Trigger spray can be set to jet or spray	<ul style="list-style-type: none"> • Large droplet size reduces the risk of inhalation and provides good surface coverage
Double bagged packed in linear tear packaging	<ul style="list-style-type: none"> • Each bag is easy to open even when wearing gloves • Facilitates transfer disinfection into cleanroom
Not classed as corrosive	<ul style="list-style-type: none"> • Can be used safely in all areas of the cleanroom

For more information or to request a sample,
e-mail infoeu@contec.com or phone + 33 (0) 2 97 43 76 98.

Copyright© 2020 Contec, Inc. All rights reserved. INT111 290220

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

Specification

- 2000ppm hypochlorous acid in purified water
- Ready-to-use with 2 year shelf life
- 0.2 micron filtered under Grade A air in a Grade B/C cleanroom
- Double bagged in polyethylene linear tear bags
- Trigger spray with jet or spray option
- Relevant certificates of analysis and sterility provided with every batch.
- Full supporting documentation information is available on request.

Efficacy

Test	Description	Log Reduction	Time	Test	Description	Log Reduction	Time
EN1276	<i>E. hirae</i>	> log 7	1 min	EN13704	<i>C. difficile</i>	> log 6	1 min
EN1276	<i>S. aureus</i>	> log 7	1 min	EN13704	<i>B. subtilis</i>	> log 6	1 min
EN1276	<i>P. aeruginosa</i>	> log 7	1 min	EN13704	<i>B. pumilis</i>	> 3.88	1 min
EN1276	<i>E. coli</i>	> log 7	1 min	EN13704	<i>B. cereus</i>	> 3.24	1 min
EN14476	Poliovirus	4.33	30 secs	EN13704	<i>P. gluconolyticus*</i>	> 3.12	1 min
EN14476	Adenovirus	4.67	30 secs	EN13697	<i>B. subtilis</i>	> log 6	1 min
EN1640	<i>A. niger</i>	> log 4	1 min	EN16615	<i>A. brasiliensis</i>	6.13 test area 1	1 min
EN1650	<i>C. albicans</i>	> log 4	1 min	EN16615	<i>A. subtilis</i>	7.11 test area 1	1 min

*House Isolate

Product Application

Ready-to-use so can be applied directly to the surface ensuring even coverage. Use an appropriate cleanroom wipe or mop and a recommended wiping technique for optimum contamination control. Leave for required contact time and wipe to dry. Contec ProChlor is low residue but where residues are an issue they are easily removed with either water or alcohol.

Part No.	Description	Packaging
SBC0502PC	Contec <i>Sterile</i> ProChlor 0.5L Capped	24 x 0.5L
SBT102PC	Contec <i>Sterile</i> ProChlor 1L Trigger Spray	6 x 1L
SBC102PC	Contec <i>Sterile</i> ProChlor 1L Capped	6 x 1L
SBC502PC	Contec <i>Sterile</i> ProChlor 5L Capped	2 x 5L
FBC0502PC	Contec ProChlor 0.5L Capped	24 x 0.5L
FBT102PC	Contec ProChlor 1L Trigger Spray	6 x 1L
FBC102PC	Contec ProChlor 1L Capped	6 x 1L
FBC502PC	Contec ProChlor 5L Capped	2 x 5L

Use biocides safely. Always read the label and product information before use.

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MAKE YOURSELF A CUPPA

2 MINS

ACHIEVE 100%
SPORICIDAL ELIMINATION 1 MIN

CONTEC PROCHLOR COMPLETE SPORICIDAL ELIMINATION IN LESS THAN 1 MINUTE*

Contec ProChlor is a new innovative sporicidal biocide with a proven log 6 kill in under 60 secs (*EN13697). Filtered and double bagged, Contec ProChlor is specifically designed for use in pharmaceutical cleanrooms. Available in a 1L trigger spray and 5L capped container.

Probably the fastest acting, cleanroom sporicide available.

Contact Contec at
infoeu@contecinc.com or by
calling 0845 652 2582
to request a sample.

www.contecinc.com

CONTEC PROCHLOR

Contec ProChlor is specifically designed for cleanroom environments. Manufactured in a controlled environment and filtered to 5 µm Contec ProChlor is especially suitable for Grade C and D cleanrooms where a fast acting and effective sporicide is needed. With no hazard classification Contec Protec has good operator acceptability and material compatibility.

WHAT IS IT ?

Contec ProChlor is a ready to use sporicidal biocide based on stabilised hypochlorous acid. Using a patented formulation Contec ProChlor provides a complete sporicidal surface kill in under 60 secs*.

WHY IS IT DIFFERENT ?

Hypochlorous acid is an equilibrium product formed when chlorine is dissolved in water. Historically, due to its very short shelf life, hypochlorous acid could only be created at point of use through the electrolysis of a salt solution. Not anymore ! Contec ProChlor is a patented stabilised version of hypochlorous acid which has a shelf life of 15 months.

WHY IS IT SO EFFECTIVE ?

It is generally accepted that hypochlorous acid is 100 times more efficacious than hypochlorite.¹ The stabilised pH of Contec ProChlor ensures that nearly all the available chlorine is in the format of hypochlorous acid creating a very fast acting and powerful sporicide which is uniquely not hazardous to people or equipment. In contrast, the available chlorine in hypochlorite solutions which is at a very high pH, is nearly all in the form of hypochlorite ions which have less sporicidal efficacy.

WHAT DO YOU GET ?

A sporicidal disinfectant which is:-

Fast acting – sporicidal on a surface in under 1 min


Cleanroom compatible – filtered and double bagged in a controlled environment

Fully validated – provided with complete testing package to EN norms


Safe to use – has no associated hazard, is not corrosive and requires no special disposal

Easy to use – supplied ready-to-use in a 1L trigger spray or 5L capped container

COMPLETE SPORICIDAL ELIMINATION

60 SECONDS LATER

For further information or for a sample of probably the fast acting, cleanroom sporicide available, contact Contec at infoeu@contecinc.com or by calling **0845 652 2582**.

* Tested to EN13704 and EN13697 for spores / ¹ White's Handbook of chlorination and alternative disinfectants / Use biocides safely. Always read the label and product information before use.