CELOS[®]

a pure¹¹ brand

CELOS GLOVE 2.30

English
German
Dutch
Spanish
Portuguese
French
Danish
Polish
Czech
Hungarian
Macedonian
Slovakian

Instruction manual

EN Nitrile Glove CELOS Glove 2.30 Category III

Sterile Nitrile Glove for Cleanroom Application

NOTE: THIS MANUAL SHOULD BE PROVIDED TO USERS OF THESE GLOVES FOR SAFETY REASONS. This product is certified as a Category III PPE product in accordance with regulation (EU) 2016/425 EEC and have been shown to comply with this Regulation through the Harmonised European Standard EN ISO 21420:2020, EN ISO 374-5:2016 and EN ISO 374-1:2016+A1:2018. SATRA Technology Europe Ltd · Bracetown Business Park Clonee · IRELAND D15 YN2P · Identification number 2777

Nitrile · Blue · Hand-specific · Powder-free · Latex-free · Acceptable Quality Level (AQL): 0.65 · Textured fingertips · Accelerator-free · Low particulate and ion residue levels · Washed with 0.2 µm filtered DI water · Triple packaging

USE: These gloves are designed to protect the user's hands from low chemical risks and microbial contamination. Gloves are part of protective clothing and protect the user from potentially toxic substances or organisms. However, always remember that no item of PPE can provide full protection and care must always be taken while carrying out risk related activities. On the other hand, they also serve to protect products and samples from contamination of any kind. These gloves are meant for single use only. Do not reuse the gloves. These gloves provide no protection against mechanical risks. The testing and labeling of the gloves comply with the PPE Regulation 2016/425, which replaces 89/686/EEC, and the applicable harmonized European standards. Please ensure that the gloves are used exclusively for their intended applications.

NOTE: This information does not indicate the actual duration of protection in the workplace or the differentiation between mixed and pure chemicals. Resistance to chemicals has been assessed under laboratory conditions on samples taken from the inside of the hand only and refers only to the chemicals tested. Verification that the gloves are suitable for the intended use is recommended, as workplace conditions may be different from those of the type test depending on temperature, abrasion and degradation. If protective gloves have already been used, they may offer less resistance to hazardous chemicals due to changes in their physical attributes. Degradation caused by contact with chemicals, movements, friction, can significantly reduce the actual application time. With aggressive chemicals, degradation may be the most important factor to consider when selecting chemical resistant gloves. According to standard EN ISO 374-5:2016: The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

STORAGE RECOMMENDATION: Do not store cleanroom nitrile gloves in locations where the temperature may exceed 104 °F (40 °C). Opened cleanroom nitrile gloves must be protected from direct sunlight or prolonged fluorescent lighting to prevent discoloration. Incorrect storage of cleanroom nitrile gloves will result in a reduced shelf life and will adversely affect their properties. Product shelf life under optimal conditions are 3 years.

ATTENTION: Primary material nitrile gloves do not contain natural rubber latex, ensuring a low allergy risk. Do not use in case of intolerance to nitrile. If an allergy was to occur, seek medical assistance immediately. Never use the product in applications which requires mechanical and/or thermal protection. Inspect gloves for damage before wearing. If there are holes, tears or re-wearing, the protective performance can no longer be guaranteed. Do not use if damaged.

DISPOSAL INFORMATIONS: The product must be disposed of in accordance with the rules and regulations in force in your country. Gloves which have been contaminated by chemical substances must be disposed of in accordance with the regulations for the chemicals concerned.

Ongoing testing for Module C2 and issuance of EU Type Examination certificate are responsible by SATRA Technology Europe Ltd, Bracetown Business Park Clonee, D15 YN2P, Ireland

You can read the conformity assessment and the issuance of the EU declaration of conformity here: <https://www.pure11.de/doc>

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DONNING AND DOFFING INFORMATION:

Donning
<ol style="list-style-type: none"> 1. Inspect the glove 2. Align the glove with your hand 3. Insert your hand 4. Check the fit

Doffing
<ol style="list-style-type: none"> 1. Grasp the outside edge of the glove 2. Peel and hold the glove 3. Slide finger into glove 4. Peel off second glove and throw away

FITTING AND SIZING: Protection is limited to the hand only. Always wear gloves that fit well and trim nails before wear. Only wear products of a suitable size. Products which are either too loose or too tight will restrict movement and will not provide the optimum level of protection. The size of these products is marked on a label attached at the packaging and carton and they are available in the following sizes: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 10.0.

Product Code	Size	Minimum Glove Length
1130311BL6	6.0	300 mm
1130311BL6_5	6.5	300 mm
1130311BL7	7.0	300 mm
1130311BL7_5	7.5	300 mm
1130311BL8	8.0	300 mm
1130311BL8_5	8.5	300 mm
1130311BL9	9.0	300 mm
1130311BL10	10.0	300 mm

EN ISO 374-5:2016	
Resistance to bacteria and fungi	Pass
Resistance to virus	Pass

No	Chemical Testing	Cas Number	Chemical Code
1	Hydrogen peroxide	7722-84-1	P
2	Formaldehyde	50-00-0	T
3	Sodium hydroxide	1310-73-2	K

Permeation Test EN ISO 374-1:2016+A1:2018			Degradation Test EN ISO 374-4:2019
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Hydrogen peroxide	> 46 - 60	Level 2	39.9

Permeation Test EN ISO 374-1:2016+A1:2018			Degradation Test EN ISO 374-4:2019
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Formaldehyde	> 480	Level 6	29.5

Permeation Test EN ISO 374-1:2016+A1:2018			Degradation Test EN ISO 374-4:2019
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Sodium hydroxide	> 480	Level 6	10.7

The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. EN ISO 374-4:2019: Degradation levels indicate the change in puncture resistance after exposure to the challenge chemical.

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SATRA Technology Europe Ltd · Bracetown Business Park Clonee · IRLAND D15 YN2P · Identifikačné číslo 2777
Importér: pure11 GmbH · Bavariafilmpatz 7 · 82031 Gruenwald · Germany