Nitril-Handschuhe Kimtech S. Sterling XTRA

pure¹¹-Nr.: 05339, Hersteller: Kimberly-Clark

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Chemical Permeation Table (printed 2021-02-16)
Standard [EN16523-1]

Glove	Concentration	CAS Number	Kimtech™ Sterling™ Xtra Nitrile Gloves	
			EN 16523-1 Permeation	EN 374-4 Degradation
Acetic Acid	99%	64-19-7	2,6	95,0%
Ammonium Hydroxide	25%	1336-21-6	5,7	75,0%
Carbon Disulphide	99%	75-15-0	<1	67,0%
Citric Acid	30%	77-92-9	>480	13,0%
Cyclohexane	99%	100-82-7	31	61,0%
Dichloromethane	100%	75-09-2	0	100,0%
Dimethyl Sulphoxide	99,50%	67-68-5	5,2	64,0%

Glove	Concentration	CAS Number	Kimtech™ Sterling™ Xtra Nitrile Gloves	
			EN 16523-1 Permeation	EN 374-4 Degradation
Ethanol	70%	64-17-5	21	65,0%
Ethidium Bromide	1%	1239-45-8	>480	5,9%
Ethyl Acetate	99%	141-78-6	<1	87,0%
Formaldehyde	37%	50-00-0	180	32,0%
Hydrochloric Acid	30%	7647-01-0	180	72,0%
Hydrogen Peroxide	30%	7722-84-1	>480	34,0%
Isopropanol	70%	67-63-0	33	67,0%
Isopropanol	99%	67-63-0	20	65,0%
Nitric Acid	65%	7697-37-2	4,0	93,0%
Perchloric Acid	70%	7601-90-3	300	0,9%
Sodium Hydroxide	40%	1310-73-2	>480	2,6%

Glove	Concentration	CAS Number	Kimtech™ Sterling™ Xtra Nitrile Gloves	
			EN 16523-1 Permeation	EN 374-4 Degradation
Sodium Hydroxide	50%	1310-73-2	>480	-8,6%
Sodium Hypochlorite	14%	7681-52-9	>480	13,0%
Tetrahydrofuran	100%	109-99-9	0	78,3%
Toluene	99%	108-88-3	<1	86,0%
n-Heptane	99%	142-82-5	20	55,0%

Disclaimer: All data provided is based on results of tests performed in accordance with the relevant test standard (Chemical Permeation: EN16523-1 or EN374-3; Degradation: EN374-4), by an independent laboratory which has approval from a notified body under the CE Regulation (Or Directive) for Personal Protective Equipment. These tests may not adequately replicate any specific conditions of use, and because KIMBERLY-CLARK PROFESSIONAL™ has no detailed knowledge or control over the conditions of end use, any of the data provided must be considered on an advisory basis only, and KIMBERLY-CLARK PROFESSIONAL* must decline any liability.


Ref: Cytotoxic Drug Permeation testing of Kimtech™ gloves

Dear Valued Customer,

Thank you for your enquiry about our gloves and their use for protection against chemical splash hazards by cytotoxic drugs.

We must caution that the selection of the most appropriate glove for a specific task should be carried out by a trained safety professional following a full risk assessment. It is Kimberly-Clark's intent to provide data which can enable trained professionals to make an informed choice.

ASTM D 6978-05 (Standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs) was utilised to complete the chemotherapy drug performance assessment on the above glove products. Testing was performed on the cuff area of the gloves, under conditions of continuous contact. UV/VIS Spectrometry was used to measure the absorbance of the challenge chemicals through the specimens, into the collection medium.

Kimtech™ Sterling™ Nitrile Xtra gloves

Chemotherapy Drug	Concentration (mg/ml)	Breakthrough Time (mins)
Carmustine	3.3	20.4
Cisplatin	1.0	No breakthrough up to 240 mins
Cyclophosphamide (Cytoxan)	20.0	No breakthrough up to 240 mins
Dacarbazine (DTIC)	10.0	No breakthrough up to 240 mins
Doxorubicin Hydrochloride	2.0	No breakthrough up to 240 mins
Etoposide (Toposar)	20.0	No breakthrough up to 240 mins
Fluorouracil	50.0	No breakthrough up to 240 mins
Ifosfamide	50.0	No breakthrough up to 240 mins
Mitoxantrone	2.0	No breakthrough up to 240 mins
Paclitaxel (Taxol)	6.0	No breakthrough up to 240 mins
Thiotepa	10.0	110
Vincristine Sulphate	1.0	No breakthrough up to 240 mins

Tests were carried out by an independent laboratory, under laboratory test conditions. These tests may not adequately replicate any specific condition of use. As KIMBERLY-CLARK* has no detailed knowledge or control over the conditions of end use, this data must be considered advisory only and KIMBERLY-CLARK* must decline any liability.

Kimberly Clark Professional
40, London Road
Reigate
Surrey
RH2 9QP
UK

Declaration of Compliance	
Manufacturer	Kimberly Clark Professional 40, London Road, Reigate, Surrey, RH2 9QP, UK
Product covered by this declaration	Kimtech Science* Sterling Nitrile Xtra Gloves
Product Codes	98341-98345
Kimberly Clark Global product Safety Clearance number	10000005181
Material Reference	N/A
Date of the Declaration	08/06/2018
Revision 00 Date	08/06/2018
Declaration Expiry date	08/06/2020
Declaration of Compliance with	
This product complies with	(EC) No. 1935/2004 (as amended)
This product complies with	(EC) No. 2023/2006 (as amended)
This product complies with	German Recommendation, BfR No. XXI (as amended)
This product complies with	(EC) No. 1907/2006 (as amended)
<p>Based upon the preceding review this article intended to come into contact with food has been formulated only with substances and additives that are authorised under German Recommendation BfR No. XXI (with reference to Regulation (EC) No. 10/2011 (as amended)).</p> <p>This article intended to come into contact with food does not contain any substance >0.1% from the list of substances of very high concern (SVHC) in the REACH Regulation (EC) No. 1907/2006.</p> <p>This article is manufactured according to regulation (EC) No. 2023/2006 (as amended) on good manufacturing practice.</p> <p>All ingredients used in this article do not exceed any legal migration levels based on intended use of the product. Raw materials used in production of the gloves are specified safe for food contact and are procured from an approved Kimberly Clark supplier.</p> <p>Kimberly-Clark has an internal product traceability system in accordance with our Quality Management System.</p>	
Article Material	
Nitrile rubber (synthetic rubber).	

Information about the compliance of substances used that are subject to any restriction or specification		
Compliance with overall migration limit	Overall migration is below 10 mg/dm ² under standard testing conditions laid down under German Recommendation, BfR No. XXI.	
Compliance with migration behavior	Migration levels are not exceeding the limitation laid down in (EC) No. 10/2011 (as amended) health – related evaluation of materials and objects for the contact with Foodstuffs in the frame of the Foodstuffs and Animal Feed Code, 13 th memorandum, Bundesgesundheitsblatt 5, 403 (1962), including the 217 th memorandum, Bundesgesundheitsblatt 57, 1350-1351 (2014), state of 1 October 2014.	
Specific migration limits (SMLs)		
The determination of Specific migration was performed on the finished article, according to the test result. No significant migration could be identified, the product not exceeding SMLs.		
Individual Substances	Specific migration Limits (SMLs)	Test results
Acrylonitrile	0.01 mg/kg	<0.01 mg/kg
Information about the compliance of substances subject to purity criteria		
There are no known substances subject to purity criteria.		
There are no known substances subject to restrictions apart from the Specific migration limits (SMLs).		
Conditions of use		
<ul style="list-style-type: none"> • Type(s) of food with which it is intended to be put in contact <ul style="list-style-type: none"> ○ All Foods. • Time and temperature and storage while in contact with the food <ul style="list-style-type: none"> ○ Maximum temperature of 40°C (104°F) up to 2 minutes. • The Ratio of food contact surface area to volume used to establish the compliance of the article <ul style="list-style-type: none"> ○ The compliance testing was done under conditions set out in German Recommendation, BfR No. XXI (as amended) using a surface to volume (s/v) contact ratio of 8.4 dm²/kg. 		
Functional Barrier		
There is no functional barrier present.		

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Erteilt an Kimberly - Clark Europe Ltd, aufgrund der technischen Unterlagen durch die benannte Stelle:
Kimberly-Clark Europe Ltd.
Kimberly-Clark Europe Limited
TÜV SÜD Product Service GmbH Zertifizierstellen (0123)