Nitril-Handschuhe Kimtech P. Sterling G3 - AUSLAUFARTIKEL

pure¹¹-Nr.: 05026, Hersteller: Kimberly-Clark

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Chemical Permeation Table (printed 2021-02-15)
Standard [EN16523-1]

| Glove | Concentration | CAS Number | Kimtech™ G3 Sterling™ Nitrile Gloves | |
|---------------------|---------------|------------|--------------------------------------|----------------------|
| | | | EN 16523-1 Permeation | EN 374-4 Degradation |
| Acetic Acid | 99% | 64-19-7 | 2,6 | 95,0% |
| Ammonium Hydroxide | 25% | 1336-21-6 | 5,7 | 75,0% |
| Carbon Disulphide | 99% | 75-15-0 | <1 | 67,0% |
| Citric Acid | 30% | 77-92-9 | >480 | 13,0% |
| Cyclohexane | 99% | 100-82-7 | 31 | 61,0% |
| Dimethyl Sulphoxide | 99,50% | 67-68-5 | 5,2 | 64,0% |
| Ethanol | 70% | 64-17-5 | 21 | 65,0% |

| Glove | Concentration | CAS Number | Kimtech™ G3 Sterling™ Nitrile Gloves | |
|---------------------|---------------|------------|--------------------------------------|----------------------|
| | | | EN 16523-1 Permeation | EN 374-4 Degradation |
| Ethidium Bromide | 1% | 1239-45-8 | >480 | 26,0% |
| Ethyl Acetate | 99% | 141-78-6 | <1 | 87,0% |
| Formaldehyde | 37% | 50-00-0 | 180 | 32,0% |
| Hydrochloric Acid | 30% | 7647-01-0 | 420 | 18,0% |
| Hydrogen Peroxide | 30% | 7722-84-1 | >480 | 34,0% |
| Isopropanol | 70% | 67-63-0 | 33 | 67,0% |
| Isopropanol | 99% | 67-63-0 | 20 | 65,0% |
| Nitric Acid | 65% | 7697-37-2 | 4,0 | 93,0% |
| Perchloric Acid | 70% | 7601-90-3 | 300 | 0,9% |
| Sodium Hydroxide | 40% | 1310-73-2 | >480 | -3,0% |
| Sodium Hypochlorite | 14% | 7681-52-9 | >480 | 13,0% |

| Glove | Concentration | CAS Number | Kimtech™ G3 Sterling™ Nitrile Gloves | |
|----------------|---------------|------------|--------------------------------------|----------------------|
| | | | EN 16523-1 Permeation | EN 374-4 Degradation |
| Sulphuric Acid | 50% | 7664-93-9 | >480 | 7,8% |
| Toluene | 99% | 108-88-3 | <1 | 86,0% |
| n-Heptane | 99% | 142-82-5 | 20 | 55,0% |

Disclaimer: All data provided is based on results of tests performed in accordance with the relevant test standard (Chemical Permeation: EN16523-1 or EN374-3; Degradation: EN374-4), by an independent laboratory which has approval from a notified body under the CE Regulation (Or Directive) for Personal Protective Equipment. These tests may not adequately replicate any specific conditions of use, and because KIMBERLY-CLARK PROFESSIONAL™ has no detailed knowledge or control over the conditions of end use, any of the data provided must be considered on an advisory basis only, and KIMBERLY-CLARK PROFESSIONAL* must decline any liability.

Ref: Cytotoxic Drug Permeation testing of Kimtech™ gloves


Dear Valued Customer,

Thank you for your enquiry about our gloves and their use for protection against chemical splash hazards by cytotoxic drugs.

We must caution that the selection of the most appropriate glove for a specific task should be carried out by a trained safety professional following a full risk assessment. It is Kimberly-Clark's intent to provide data which can enable trained professionals to make an informed choice.

ASTM D 6978-05 (Standard practice for assessment of resistance of medical gloves to permeation by chemotherapy drugs) was utilised to complete the chemotherapy drug performance assessment on the above glove products. Testing was performed on the cuff area of the gloves, under conditions of continuous contact. UV/VIS Spectrometry was used to measure the absorbance of the challenge chemicals through the specimens, into the collection medium.

K-C Confidential. The information contained in this document is intended for the use of the addressee(s) only and may contain privileged, confidential, or proprietary information that should not be distributed publically without Kimberly-Clark's written consent.

Kimtech™G3 Sterling™ Nitrile gloves

| Chemotherapy Drug | Concentration (mg/ml) | Breakthrough Time (mins) |
|----------------------------|-----------------------|--------------------------------|
| Carmustine | 3.3 | 20.4 |
| Cisplatin | 1.0 | No breakthrough up to 240 mins |
| Cyclophosphamide (Cytoxan) | 20.0 | No breakthrough up to 240 mins |
| Dacarbazine (DTIC) | 10.0 | No breakthrough up to 240 mins |
| Doxorubicin Hydrochloride | 2.0 | No breakthrough up to 240 mins |
| Etoposide (Toposar) | 20.0 | No breakthrough up to 240 mins |
| Fluorouracil | 50.0 | No breakthrough up to 240 mins |
| Ifosfamide | 50.0 | No breakthrough up to 240 mins |
| Mitoxantrone | 2.0 | No breakthrough up to 240 mins |
| Paclitaxel (Taxol) | 6.0 | No breakthrough up to 240 mins |
| Thiotepa | 10.0 | 110 |
| Vincristine Sulphate | 1.0 | No breakthrough up to 240 mins |

Tests were carried out by an independent laboratory, under laboratory test conditions. These tests may not adequately replicate any specific condition of use. As KIMBERLY-CLARK* has no detailed knowledge or control over the conditions of end use, this data must be considered advisory only and KIMBERLY-CLARK* must decline any liability.

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¹ Nitril ist ein synthetisches Material, das viele Eigenschaften mit Naturkautschuk-Latex gemein hat, sich aber von diesem durch mehrere signifikante Vorteile unterscheidet: hoher Tragekomfort, gute Stichfestigkeit, hohe Abriebfestigkeit ohne Beeinträchtigung der Tastempfindlichkeit oder der elektrostatisch dissipativen Eigenschaften. ² AQL-Wert gemäß Definition nach ISO 2859-1 für Attributprüfung anhand von Proben. ®/™ Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. © KCWW. The colour Gray and Sterling™ are trademarks of Owens & Minor, Inc. or its affiliates and used under license. Publication code: ID4483.04 DE 07.20
Erteilt an Kimberly - Clark Europe Ltd, aufgrund der technischen Unterlagen durch die benannte Stelle: