

## Chemoprotect Arbeitsunterlage

pure<sup>11</sup>-Nr.: 06407, Hersteller: Distributor pure<sup>11</sup>



### Zusammenfassung

- Neue pure<sup>11</sup>-Artikelnummer (ab 01.07.2023): 1106407
- Farbe: Blau
- Einsatz: Bei der Zubereitung von Zytostatika
- Bestehend aus 3 Lagen
- Sterilisiert durch Ethylenoxid
- Einzeln-verpackt (doppelt)

### Empfohlene Reinraumklassen

ISO	3	4	5	6	7	8	9
GMP			A/B		C	D	

### Produktvarianten

pure<sup>11</sup>-Nr.: 06407


Maße: 41 x 54 cm / VE: 50 Stück

Quelle: <https://www.pure11.de/chemoprotect-arbeitsunterlage-0>



## Chemoprotect® Prep mat

### Chemoprotect® Prep mat - sterile



1. The Chemoprotect® Prep mat has been especially developed for the preparation of cytotoxic drugs. The three-layer structure of the material allows safe working and offers a workplace that is non-slipping and free from fluffs. The Chemoprotect® Prep mats are very well suited for isolators respectively safety work benches. The liquid absorption of Chemoprotect® Prep mat minimises the risk of contaminations during regular handling of cytotoxic drugs as well as in case of spillage. Furthermore, it safeguards a subsequent disposal.

#### Product Specification

##### Material:

- Lint-free fluid-permeable upper layer surface (Polyethylene)
- Highly absorbent fluid-retaining middle layer (non-woven material), bleached virgin wool pulp
- Non-slip lower layer, impermeable to liquids (Polyethylene)

##### Features:

- Liquid-permeable surface, absorbent, liquid-receptive middle layer, liquid impermeable lower layer and non-fuzzing

##### Absorption of liquids:

- $\leq 950 \text{ ml/m}^2$  (absolute)
- $\leq 100 \text{ ml/Chemoprotect}^{\circledR}$  Prep mat, small
- $\leq 200 \text{ ml/Chemoprotect}^{\circledR}$  Prep mat, large



#### Disclaimer:

The above is indicated only for information. This document is not a specification. The values given are subject to the usual fluctuations. Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

Subject to alterations without prior notice!

REF	Product description	Packaging	GTIN	Filling/Storage volume
79.0301	Chemoprotect® Prep mat small (27 x 42 cm) - sterilised	1 50 300 1800	04055698013463 04055698013456 04055698013470 -	



## Chemoprotect® Prep mat

### Chemoprotect® Prep mat - sterile


REF	Product description	Packaging	GTIN	Filling/Storage volume
79.0303	Chemoprotect® Prep mat large (41 x 54 cm) - sterilised	1 50 300 1800	04055698013487 04055698013494 04055698013500 -	

## Chemoprotect<sup>®</sup> Prep mat, large - sterile

REF 79.0303

### 10 Prep mat

- Film: polyethylene film treated for anti-static properties
- Paper: airlaid nonwoven absorbent material made from bleached virgin wood pulp
- Adhesive: hot melt
- Netting: polyethylene



### Disclaimer:

Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

**Subject to alterations without prior notice!**

### Chemoprotect® Prep mat


The Chemoprotect® Prep mat has been especially developed for the preparation of cytotoxic drugs. The three-layer structure of the material allows safe working and offers a workplace that is non-slipping and free from fluffs. The Chemoprotect® Prep mats are very well suited for isolators respectively safety work benches. The liquid absorption of Chemoprotect® Prep mat minimises the risk of contaminations during regular handling of cytotoxic drugs as well as in case of spillage. Furthermore, it safeguards a subsequent disposal.



### Product Specification

#### Material:

- (1) Lint-free fluid-permeable upper layer surface (Polyethylene)
- (2) Highly absorbent fluid-retaining middle layer (non-woven material, bleached virgin wood pulp)
- (3) Non-slip lower layer, impermeable to liquids (Polyethylene)



#### Features:

Liquid-permeable surface, absorbent, liquid-receptive middle layer, liquid impermeable middle layer and non-fuzzing

#### Absorption of liquids:

≤ 950 ml/m<sup>2</sup> (absolute)  
≤ 100 ml/Chemoprotect® Prep mat, small  
≤ 200 ml/Chemoprotect® Prep mat, large



### Product range - non-sterile

REF	Description	Size	Packaging units
79.0300	Chemoprotect® Prep mat, small	27 x 42 cm	100 pcs. /dispenser box 1000 pcs. /shipping box
79.0302	Chemoprotect® Prep mat, large	41 x 54 cm	50 pcs. /dispenser box 500 pcs. /shipping box

### Product range - sterile

REF	Description	Size	Packaging units
79.0301	Chemoprotect® Prep mat, small	27 x 42 cm	50 pcs. /dispenser box 300 pcs. /shipping box
79.0303	Chemoprotect® Prep mat, large	41 x 54 cm	50 pcs. /dispenser box 300 pcs. /shipping box

#### Please note:

The above mentioned products do not comply with at least 1 of approx. 20 EC directives so that a CE mark cannot be affixed! This also applies with consumables such as foils tubes used for the purpose to weld medicine or waste or aid-consumables such as absorbent- or wiping tissues for liquids.

#### Disclaimer:

The above is indicated only for information. This document is not a specification. The values given are subject to the usual fluctuations. Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

**Subject to alterations without prior notice!**